


Ridge Expander Kit

Unit: mm, Scale 1 : 1




GREK

Ridge Spreader

Diameter	L	Art No.
Ø1.3 / Ø2.5	32	GRS 13 25
Ø1.8 / Ø3.0	32	GRS 18 30
Ø2.3 / Ø3.5	32	GRS 23 35


※ Note: condensing drill 20~60 rpm,
30~45N-cm with irrigation



Harvest Drill


Diameter	L	Art No.
Ø3.0	29	GHD 30 29

※ Note: Harvest Drill speed 20~60rpm,
30~45 N-cm without irrigation



Disc Saw

Diameter	L	Art No.
Ø8.0	25	GDS 80 25



Bone Chisel

GBC 18 45 13

