

bright 3-Layer

Minimal Coloring & Staining

3 Layer

Shade

A1, A2, A3, A3.5

Thickness

14T, 16T, 18T, 22T, 30T



Incisal Layer (4.3mm)


4Y/5Y Zirconia (650MPa, 49%)

Body Layer (4.3mm)

4Y/5Y Zirconia

Cervical Layer (remaining thickness)

4Y Zirconia (950MPa, 46%)



Recommended

Anterior
Crown
Single Crown



Recommended

Posterior
Crown
Single Crown



Recommended

Anterior
Bridge
Bridge (3-unit-Bridge)




Recommended

Posterior
Bridge
Bridge (3-unit-Bridge)




Bridge (≥4-unit-Bridge)



Positioning



Coloring



49% Transparency

Final prothesis